


FDA'S Interlaboratory Computational Study Of Flow In A Nozzle Model: Using Experimental Results To Refine Simulations

Sandy F.C. Stewart¹, Prasanna Hariharan¹, Greg W. Burgreen²,
Eric G. Paterson³, Varun Reddy³, Steven W. Day⁴, Matthew Giarra⁴,
Keefe B. Manning³, Steven Deutsch³, Matthew R. Myers¹,
Michael R. Berman¹, Richard A. Malinauskas¹


¹Food & Drug Administration, ²Mississippi State University,
³Pennsylvania State University, ⁴Rochester Institute of Technology

Computational Round Robin



- 5 flow rates with throat Reynolds number (Re_t) = 500, 2000, 3500, 5000, and 6500
- 28 groups submitted simulation results (including 3 from organizing committee)
- Experimental validation performed via particle image velocimetry (PIV) at 3 labs

Experimental Validation




- Blood analog fluid: 50% NaI, 20% Glycerin, 30% water
- Refractive index matches acrylic model
- PIV performed at same Reynolds numbers as CFD
- PIV results scaled to match CFD fluid properties


Validating CFD to Experiment

- Asked participants to redo simulations & match results to experimental data
- Sudden Expansion
- $Re = 500$ and 3500
- 10 outside groups provided redo simulations
- 3 in organizing committee also redid simulations


$Re_t = 500$, Centerline Velocity




$Re_t = 500$, Velocity Profile (0.06 m)



$Re_t = 3500$, Centerline Velocity




$Re_t = 3500$, Velocity Profile (0.016 m)




Simulation Redo – Internal Study

$Re_t = 3500$, Centerline Velocity




Simulation Redo – Internal Study

$Re_t = 3500$, Velocity Profile (0.016 m)




Mesh Used in Internal Study Redo

$Re_t = 3500$ (k- ω , fine mesh, no wall functions)



Convergence in Internal Study Redo

$Re_t = 3500$ (k- ω , fine mesh no wall functions)




A Few Best Practices

- Choose appropriate physics
 - Choose the right turbulence model if needed
- Choose appropriate boundary conditions and inlet/outlet lengths
- Ensure conservation of mass
- Grid refinement studies are necessary but must be combined with other good practices
- In the transitional regime, experimental results may not be unique

A Few Best Practices

- Transitional cases may require time-dependent or LES model to fully capture details
- Good methods from **this** study may not be transferable to **other** devices/flow regimes
- Good convergence criteria is essential
- VALIDATION
- VALIDATION
- VALIDATION

Future of FDA Project



- Blood damage testing under way at 3 labs for comparison
- VAD prototype is being tested for future interlaboratory CFD/PIV study
- FDA to work with standards group (AAMI, ASME)

Acknowledgments

We sincerely thank our 49 study participants, from 6 countries

Yared Alemu	Marcus Hormes	T. Schauer
K. Amano	Marc Horner	E. Sirois
J. Ashton	Xueying Huang	Ken Solen
Mehdi Behbahani	S. Krishnan	Ulrich Steinseifer
Catrin. Bludszuweit-Philipp	S. Kühne	Wei Sun
Danny Bluestein	Patricia Lawford	S. Takahashi
Parnian Boloori-Zadeh	Ming Li	Masaaki Tamagawa
Alistair Brown	Z. Li	K. Tan
Scott Corbett	William Louisos	Dalin Tang
Julien de Charentenay	Takehisa Mori	Ertan Taskin
D. de Zelicourt	Alex Medvedev	Z. Teng
Wade Dummer	Joerg Mueller	R. Toyao
K. Fraser	Sei Murakami	Dan White
Leonid Goubergrits	Mariagrazia Pilotelli	Z. Jon Wu
Kurt Graichen	M. Polster	M. Xenos
Linden Heflin	N. Reuel	Ajit Yoganathan
Darren Hitt		